

Table of Contents

Introduction	2
Figure 1	2
Thoughts or questions.	2
Protocol 1	2
Reagents	3
File Formatting Protocol	3
September 2018	4
Bibliography	5

Introduction

Thursday, September 20, 2018

To Do:

- ~~1. Explain lab notebook formatting~~
2. Media prep

Results and Data:

For each day, copy the previous entry headers and update the date. Save the same ongoing copy of the lab notebook to the lab google drive ONLY. Do not save separate physical copies. At the 1st of each new month, a read-only PDF file of the lab notebook must be saved and given to Dr. Ramsey electronically.

File contents converted to PDF MUST NOT BE EDITED after PDF conversion. Continue to keep using the same word file until the end of the calendar year. New lab notebook files run in 6 month intervals.

For the To do list, update this each day with new tasks, as tasks are done, use the ~~strike through~~ font on the day they are completed and leave them out of the list on the next day.




Figure 1

Figures are inserted as inline .png files when possible, .jpg, .gif, .tif otherwise are acceptable. Figure legends are always inserted in Word (right click – ‘insert caption’) and use “**Heading 2**” text to properly format in the TOC. Table legends are handled the same way. Figure images must be saved in a separate folder where the source data is saved.

Thoughts or questions. When you have a significant observation, question, confusing point or contradiction that you have identified in your data or protocol, use the “**Heading 3**” text heading on a descriptive brief title or single word heading so you can refer to it in your TOC later.

Protocol 1

1. This is the format for a protocol in your notebook.

2. The protocol title must be formatted in “**Heading 2.**”
3. Reagents which must be made for a specific protocol (buffers, solutions) should be listed in a “Reagents” Section, formatted in “**Heading 3.**”
4. The protocol must be in numerical steps.
5. Use standard notation and carefully describe units for your protocol.
6. Use ½” indent for protocol text.

Reagents

Specific buffers

For initial lab notebooks, write as much detail as possible. As time goes by you will be able to refer to written protocols by their heading and/or page number within the document. If you make any modifications to a protocol you must state how and why.

File Formatting Protocol

1. Filenames begin with your initials, an underscore, and the date, formatted as the last 2 digits of the year, the month, then the day, ex: “KMR_180920_Sample_file1v1.xls”
2. This ensures all files will be sorted by their creator and by their date. You must use this file formatting system for all data files (including photos) that will be shared with the lab.

Bibliography data will be saved as (author/date) and using Mendeley at this time with the TOC inserted by order cited at the end of the lab notebook in the Bibliography section.

Example is a recent publication (Ramsey and Dove, 2016).

October 2018

Monday, October 22, 2018

To Do:

1. Prep for LB Media

LB Broth Protocol 500 mL

1. Weigh out and add into 1L Graduated cylinder/Flask
 - 10g Agar
 - 5g NaCl
 - 5g Tryptone
 - 2.5g Yeast extract
2. Add 500mL diH₂O (Don't forget to add stir bar)
3. Stir solution in hot plate until solutes dissolve
4. Cover cylinder/flask with tinfoil and stick it with autoclave tape with date written out
5. Autoclave under Liq 30°C cycle setting
6. Pour media onto plates using good microbiological technique

Tuesday, October 23, 2018

To Do:

- ~~1. Make LB Media~~
- ~~2. Sterilize stick and culture tubes~~

Results and Data:

LB Broth Protocol 1000 mL

1. Weigh out 10g of tryptone into a 1L graduated cylinder/Flask
2. Add 10g NaCl and 5g yeast extract
3. Add ddiH₂O (type 1) 1 liter
4. Mix until dissolved add stir bar
5. Aliquot out 3x200 mL 1L flasks and 4x75 mL in 100mL bottles
6. Autoclave on 30 C liquid cycle, with water in the bin

Making Chemically Competent Cells Protocol

Inoculate overnight culture with a single colony in 3mL of LB plus appropriate antibiotic

Incubate culture in 37 °C overnight

Add 3mL of filter sterilized 1M MgCl₂ to 200 mL of LB, also add appropriate antibiotic

Use 0.5 mL of overnight culture to inoculate the 200mL of LB

Incubate culture at 37 °C with aeration until OD600 of 0.5 is reached, approximately 4 hrs for DH5 α cells.

Transfer cells to a large sterile centrifuge bottle

Pellet cells by centrifugation at 4°C

Resuspend cell pellet in 60mL of cold solution A

Incubate suspension on ice for 20 min

Pellet cells by centrifugation again at 4 and resuspend in 12mL of cold solution A containing 15% glycerol

Aliquot competent cells in 0.5-1.0 mL volumes in sterile pre-cooled 1.5mL microcentrifuge tubes

Quick freeze cells using dry ice

Smaller scale:

Add 300 μ L of filter sterilized 1M MgCl₂ to 20 mL of LB, also add appropriate antibiotic

Use 50 μ L of overnight culture to inoculate the 20mL of LB

Incubate culture at 37 °C with aeration until OD600 of 0.5 is reached, approximately 3 hrs.

Transfer cells to falcon tube or 2mL eppi tubes

Pellet cells by centrifugation at 4°C

Resuspend cell pellet in 6mL of cold solution A

Incubate suspension on ice for 20 min

Pellet cells by centrifugation again at 4 and resuspend in 1.2mL of cold solution A containing 15% glycerol

Aliquot competent cells in 500-100 μ L volumes in sterile pre-cooled 1.5mL microcentrifuge tubes

Store at -80(flash freeze optional)

Solution A:

Combine 10mL of 1M MnCl₂, 50 mL 1M CaCl₂, 200 mL 50 mM 2-morpholino-ethanesulfonic acid (MES) pH 6.3, and 740 mL ddH₂O

Solution A %15 glycerol

10mL of 1.0 MnCl₂

50mL 1.0M CaCl₂

200mL MES pH 6.3

590 mL of ddiH₂O

150mL of glycerol

Need 250mL of MES

MES= 213.25g/mol

Wednesday, October 24, 2018

To Do:

- ~~1. Streak E. coli XL1 Blue on LB plates~~

Thursday, October 25, 2018

To Do:

- ~~1. Record results of E. coli XL1 Blue streak plate~~

Results:

Isolated colonies present



Friday, October 26, 2018

To Do:

1. ~~Make 1L 50mM MES solution pH (6.3)~~
2. ~~Make 1L Solution A and 1L Solution A + 15% Glycerol~~

Monday, October 29, 2018

To Do:

1. PCR of oligos positive, negative, and experimental controls (p822 For and p823 Rev) 304 bp region amplification region
- control primers: p780 and p781 (443 bp range)

Standard PCR Protocol

1. Add ddiH₂O to 1.5mL tube
2. Add KOD buffer
3. Add KOD
4. Mix by pipetting
5. Aliquot into PCR tubes tube#1 1.5μ of p822 and 1.5μL of p823
6. Aliquot controls into tube 2 and tube 3

Total reaction volume		50		
Total number		3		

of reactions				
				Factor
Component	Stock concentration	Final concentration	1 rxn volume	4.4
ddiH ₂ O			10	44
KOD buffer	2x	1x	25	110
dNTPs	2 mM	0.4 mM	10	44
oligo F	10 uM	0.3 uM	1.5	.
oligo R	10 uM	0.3 uM	1.5	
template	100 ng/ul	2 ng/ul	1	4.4
KOD	1 U/ul	0.02 U/ul	1	4.4
		Total volume	50	220

Tuesday, October 30, 2018


To Do:

- ~~1. Run 1% Agarose gel with PCR products of 10/29/18~~
- ~~2. Visualize gel~~

Expected Size

Lane 1	Contents 2 Log NEB Ladder	Varies
Lane 2	Contents CS#1	904 bp
Lane 3	Contents CS#2	781bp
Lane 4	Contents CS#3	0 bp
Lane 5	Contents JC#1	p824/p825 824bp
Lane 6	Contents JC#2	443 bp
Lane 7	Contents JC#3	
Lane 8	Contents KR#1	
Lane 9	Contents KR#2	
Lane 10	Contents KR#3	

Results:



Wednesday, October 31, 2018

To Do:

1. ~~Make O/N of XL1-Blue~~
2. ~~Restructure Lab Notebook.~~

November 2018

Thursday, November 1, 2018

Making Chemically Competent E. Coli

1. Add 3 mL of sterile $MgCl_2$ to 200 mL LB in 1 L flask. Add antibiotic if appropriate.
2. Inoculate 200 mL LB with 0.5 mL of culture grown overnight.
3. Incubate 200 mL culture at 37°C, shaking until culture reaches an OD₆₀₀ of 0.5. ~6 hours incubation time yielded a OD₆₀₀ of .55
4. Monitor culture growth by assessing OD₆₀₀ using the spectrophotometer: At 3, 4, and 5 hours or at appropriate times between, measure OD₆₀₀
Use a cuvette with 1 mL LB as a blank
Add 1 mL culture from flask to a cuvette
Use MRamsey lab spectrophotometer set to a wavelength of 600nm
5. When cultures approach correct OD, cool down centrifuge to 4°C (Dutta lab benchtop centrifuge)
6. When culture reaches an OD₆₀₀ of approximately 0.5, transfer culture volume to sterile tubes (4x 50 mL conical) to pellet bacteria
7. Place tubes in cool centrifuge and pellet bacteria by spinning (in Dutta lab benchtop centrifuge, 15 minutes at 4000 rpm).
8. Remove tubes from centrifuge, decant the supernatant into a waste bottle, and keep cell pellets on ice
9. Add a total of 60 mL cold solution A to cell pellets from the original 200 mL culture. If you are using four 50 mL conical tubes, each tube should contain 15 mL cold solution A.

10. VERY gently resuspend cell pellet by pipetting up and down. Don't completely dispense liquid with each cycle, to prevent creating bubbles/froth. The cells should be completely homogeneous when done (no clumps or chunks)
11. Incubate resuspended cells on ice for at least 20 minutes (can stay on ice for up to 3 hours).
12. Place tubes in cool centrifuge and pellet bacteria by spinning (in Dutta lab benchtop centrifuge, 15 minutes at 4000 rpm).
13. While cells are spinning, prepare tubes for final competent cell aliquots: label and pre-chill on ice.
14. Remove tubes from centrifuge, decant the supernatant into a waste bottle, and keep cell pellets on ice
15. Add a total of 12 mL cold solution A + 15% glycerol to cell pellets from the original 200 mL culture. If you are using four 50 mL conical tubes, each tube should contain 3 mL cold solution A + 15% glycerol.
16. VERY gently resuspend cell pellet by pipetting up and down as previously.
17. Aliquot competent cells in 550 uL volumes into sterile pre-cooled 1.5 mL microcentrifuge tubes.
18. If available, freeze cells immediately upon aliquoting using using dry ice.
19. Store competent cells at -80°C.
20. At first use, test competency of cells by transforming with a known amount of supercoiled plasmid and record the transformation efficiency.

Tuesday, November 6, 2018

To Do:

1. ~~Transform chemically competent *E. coli* cells~~

Tube number	Purpose	DNA	Volume of DNA	Final volume to plate	Number of kanamycin-containing plates	Number of carbenicillin-containing plates
1	(+) control	pKL02	1 uL	20 ul, 100 ul, remaining	3	0
2	(+) control	pKL80	1 uL	20 ul, 100 ul, remaining	3	0
3	(+) control	pF	1 uL	20 ul, 100 ul, remaining	3	0
4	(-) control	None	0	20 ul, 100 ul, remaining	3	0
5	(+) control	pUC19	1 uL	20 ul, 100 ul, remaining	0	3
6	(-) control	None		20 ul, 100 ul, remaining	0	3
Total number of plates					12	6

Reaction Table

1. Check to be sure you have enough plates with appropriate antibiotic. If plates were stored at 4°C, warm at 37°C until needed. (Blue Streak=Kanamycin Black Streak= Carbenicillin)
2. Obtain DNA and thaw on ice if necessary.
3. Thaw appropriate number of competent cell tubes on ice (5 reactions per tube of competent cells)
4. Label sterile tubes as indicated in reaction table. Add indicated volume of indicated DNA on ice.

5. When competent cells are thawed (check by probing for frozen cells using a sterile pipette tip), gently pipette 100 μ L of cells into each reaction tube directly onto DNA.
6. Incubate cells on ice for 20 minutes. During incubation, find or set heat block to 42°C.
7. Place tubes with cells and DNA onto 42°C heat block for 30 seconds (heat shock step).
8. After heat shock, place tubes back on ice until next step (don't keep them here too long).
9. Using aseptic technique, add 1 mL LB (no antibiotic) to each microfuge tube.
10. Using autoclave tape, tape microfuge tubes down in shaking incubator set to 37°C.
11. Allow cells to recover for 1 hour at 37°C, shaking.
12. Using aseptic technique, plate indicated amount of cells on appropriate antibiotic plates, spreading until plates look dry. For "remaining" volume, spin tubes at max speed in benchtop centrifuge for 30 seconds (**Tube 4 & Tube 6 were not spun**). Remove 800 μ L of media. Using 200 μ L pipette, resuspend cells at bottom of tube and plate all the remaining culture. Then incubate @ 37 degrees O/N

Wednesday, November 7, 2018

To Do:

- ~~1. Count colonies on spread plates~~
- ~~2. Calculate Transformation Efficiencies~~
- ~~3. Inoculate 5mL tube with colonies of the pKL02, pKL80, and pF constructs~~

Plasmid Construct					
Volumes Aliquoted	pKL02	pKL80	pF	PUC19	None
20 μ L	7	39	7	0	0
100 μ L	80	290	122	8	0
Remaining ~800 μ L	TMTC	TMTC	TMTC	17	0
Transformation Efficiency	20 μ L=	20 μ L=	20 μ L=	20 μ L=0	20 μ L=0
	100 μ L=	100 μ L=	100 μ L=	100 μ L=1.76x10 ⁶	100 μ L=0
	Remaining=	Remaining=	Remaining=	Remaining=6.12x10 ⁶	Remaining=0

Thursday, November 8, 2018

To Do:

1. ~~Miniprep XL1 Blue constructs~~

Wednesday, November 14, 2018

Guanine Quadruplexes

To Do:

1. ~~Run LVS Genome Fasta file through G4 IPDB web tool~~
2. ~~Create Excel sheet of LVS G4 IPDB results~~

Results:

Refer to 18.11.14 LVS Genome G4 Predictor v.2 (XLSX File) Date modified Dec.4.18 in Team Google Drive - KRamsey Lab – Dan Rosario

<https://drive.google.com/open?id=1ZMhzIF0yH9DvGdUipqSWmEm-K3xd5AyJ>

Thursday, November 15, 2018

To Do:

1. ~~Create BED Files of LVS G4 IPDB files~~

Results: Time in laboratory was used on creating Bed File for Antisense file.

Refer to 181115_LVSG4_Predictor_AntiSenseStrand.bed in Team Google Drive – Kramsey Lab – Dan Rosario

https://drive.google.com/open?id=1ZkzjpjALFxR_CoJpauRn1YPJ1hSxQqn3

Week of November 18, 2018 - November 24, 2018

Thanksgiving Grace Period - NO LAB

Week of November 25, 2018 - December 1, 2018

LAB Cancellation

Reason: Dr. Ramsey's Sick Days

December 2018

Tuesday, December 4, 2018

To Do:

1. Practice and become familiar with unix/linux command terminus Task: Use BED tools on created LVS promoter, G4 Prediction antisense BED files
2. Create worklist/work document of command prompts
3. Enlist help to install linux on personal laptop

Results:

18.12.04

I tried to use the help command:

143:~ kathrynamsey\$ bedtools --help

bedtools is a powerful toolset for genome arithmetic.

Version: v2.27.1
 About: developed in the quinnlab.org and by many contributors worldwide.
 Docs: <http://bedtools.readthedocs.io/>
 Code: <https://github.com/arq5x/bedtools2>
 Mail: <https://groups.google.com/forum/#!forum/bedtools-discuss>

Usage: bedtools <subcommand> [options]

The bedtools sub-commands include:

[Genome arithmetic]

intersect	Find overlapping intervals in various ways.
window	Find overlapping intervals within a window around an interval.
closest	Find the closest, potentially non-overlapping interval.
coverage	Compute the coverage over defined intervals.
map	Apply a function to a column for each overlapping interval.
genomcov	Compute the coverage over an entire genome.
merge	Combine overlapping/nearby intervals into a single interval.
cluster	Cluster (but don't merge) overlapping/nearby intervals.
complement	Extract intervals _not_ represented by an interval file.
shift	Adjust the position of intervals.
subtract	Remove intervals based on overlaps b/w two files.
slop	Adjust the size of intervals.
flank	Create new intervals from the flanks of existing intervals.
sort	Order the intervals in a file.
random	Generate random intervals in a genome.
shuffle	Randomly redistribute intervals in a genome.
sample	Sample random records from file using reservoir sampling.
spacing	Report the gap lengths between intervals in a file.
annotate	Annotate coverage of features from multiple files.

[Multi-way file comparisons]

multiinter	Identifies common intervals among multiple interval files.
unionbedg	Combines coverage intervals from multiple BEDGRAPH files.

I tried to use the Bedtools intersect command

```
bedtools intersect -wb -a Promoters.bed -b 181115_LVSG4_Predictor_AntiSenseStrand.bed.txt >
181204_Promoter_Gpredict_antisense1.bed
```

```
bedtools intersect -wb -a Promoters.bed -b 181115_LVSG4_Predictor_AntiSenseStrand.bed.txt >
181204_Promoter_Gpredict_Sense1.bed
```

<https://drive.google.com/open?id=1axV7-DX-Szct3t59DtHNF3kDf02Dgvc5xI3wQ0Md3Vs>

Wednesday, December 5, 2018

To Do:

1. Compare Promoter Sense and Antisense G4 prediction files for sequence overlaps

Results:

18.12.05

Cat 181204_Promoter_GPredict_Sense1.bed | head

```
NC_007880 10945 10964 NC_007880 10945 10964 100
NC_007880 19346 19371 NC_007880 19346 19371 90
NC_007880 35824 35852 NC_007880 35824 35852 100
NC_007880 39598 39619 NC_007880 39598 39619 100
NC_007880 39749 39775 NC_007880 39749 39775 90
NC_007880 59174 59197 NC_007880 59174 59197 90
NC_007880 61646 61674 NC_007880 61646 61674 110
NC_007880 67877 67904 NC_007880 67877 67904 110
NC_007880 79978 80004 NC_007880 79978 80004 100
NC_007880 81311 81339 NC_007880 81311 81339 80
```

cat 181204_Promoter_GPredict_Sense1.bed | wc -l

169

cat 181204_Promoter_GPredict_antiSense1.bed | wc -l

Dan Rosario kathrynramsey\$ cat 181204_Promoter_GPredict_antiSense1.bed | wc -l

169

The 181204_Promoter_GPredict_antiSense1.bed file and the 181204_Promoter_GPredict_Sense1.bed file yielded the same output of 169 overlaps. The identical results indicate that the bed intersect command prompts were mislabeled.

```
bedtools intersect -wb -a Promoters.bed -b 181115_LVSG4_Predictor_AntiSenseStrand.bed.txt >
181204_Promoter_GPredict_antisense1.bed
```

```
bedtools intersect -wb -a Promoters.bed -b 181115_LVSG4_Predictor_AntiSenseStrand.bed.txt >
181204_Promoter_GPredict_Sense1.bed
```

Yield of same output of 169 overlaps for both strands suggest that the bed intersect command prompts were mislabeled.

Correction

```
bedtools intersect -wb -a Promoters.bed -b 181204_LVSG4_Predictor_SenseStrand.bed >
181204_Promoter_Gpredict_Sense1.bed
```

18/12/06

Command line bedtools

Output: Bedtools feature data

```
daniel@daniel-Virtualbox:~$ ls
change directory (up one directory)
```

```
daniel@daniel-Virtualbox:~$ cd ..
daniel@daniel-Virtualbox:~$ ls
```

```
daniel@daniel-Virtualbox:~$ cd ..
```

```
daniel@daniel-Virtualbox:~$ ls
```

```
daniel@daniel-Virtualbox:~$ cd ..
```

```
daniel@daniel-Virtualbox:~$ cd home/daniel/
daniel@daniel-Virtualbox:~$ ls
```

Ideally it would be beneficial to copy, paste, and or cut terminal command line inputs from the terminal itself to the main computer outside Virtualbox(Ubuntu) for the purpose of documentation. After many attempts of trying to copy/paste between the host machine and the virtualbox machine I decided to try copy pasting inside of the virtualbox

1. I logged into my google drive to access this note file using the Firefox web browser within Ubuntu
2. Copied the command lines that allowed me to change directory
3. I then returned to this note and was able to ctrl v the commands and their output. Commands and their outputs are listed below

Out of ease and convenience I can perform all computational work within the Ubuntu machine in Virtualbox and document it.

```
daniel@daniel-VirtualBox:~$ ls
181204_Promoter_Gpredict_antisense1.bed Downloads Pictures Videos
Desktop examples.desktop Public
Documents Music Templates
daniel@daniel-VirtualBox:~$ cd ..
daniel@daniel-VirtualBox:/home$ ls
daniel
daniel@daniel-VirtualBox:/home$ cd ..
daniel@daniel-VirtualBox:/$ ls
bin dev initrd.img lib64 mnt root snap sys var
boot etc initrd.img.old lost+found opt run srv tmp vmlinuz
```

```

cdrom home lib      media  proc sbin swapfile usr vmlinuz.old
daniel@daniel-VirtualBox:/$ cd ..
daniel@daniel-VirtualBox:/$ ls
bin dev initrd.img  lib64  mnt root snap  sys var
boot etc initrd.img.old lost+found opt run srv  tmp  vmlinuz
cdrom home lib      media  proc sbin swapfile usr vmlinuz.old
daniel@daniel-VirtualBox:/$ cd ..
daniel@daniel-VirtualBox:/$ cd home/
daniel@daniel-VirtualBox:/home$ cd daniel/

```

18.12.07

Changing terminal directories to “Documents”. “Documents” contains the working bed files. The transcript from the terminal below shows my attempts at changing the directories. <http://www.linuxandubuntu.com/home/10-basic-linux-commands-that-every-linux-newbies-should-remember>

The web link above is an indispensable resources for teaching basic linux commands. My rationale for using for “mkdir” which stands for make directory was ultimately not necessary since “Documents” was already a directory with the bed files and the ramseylabwork folder within it.; therefore ignore this part of the transcript altogether. cd was utilized to chnge to documents but “documents” is the incorrect directory. “cd Documents” proved to be the correct command and was verified by “ls”

```

daniel@daniel-VirtualBox:~$ mkdir ramseylabwork
daniel@daniel-VirtualBox:~$ mkdir
mkdir: missing operand
Try 'mkdir --help' for more information.
daniel@daniel-VirtualBox:~$ mkdir help
daniel@daniel-VirtualBox:~$ mkdir --help
Usage: mkdir [OPTION]... DIRECTORY...
Create the DIRECTORY(ies), if they do not already exist.

```

Mandatory arguments to long options are mandatory for short options too.

- m, --mode=MODE set file mode (as in chmod), not a=rwx - umask
- p, --parents no error if existing, make parent directories as needed
- v, --verbose print a message for each created directory
- Z set SELinux security context of each created directory to the default type
- context[=CTX] like -Z, or if CTX is specified then set the SELinux or SMACK security context to CTX
- help display this help and exit
- version output version information and exit

GNU coreutils online help: <<http://www.gnu.org/software/coreutils/>>

Full documentation at: <<http://www.gnu.org/software/coreutils/mkdir>>


```

or available locally via: info '(coreutils) mkdir invocation'
daniel@daniel-VirtualBox:~$ mkdir documents
daniel@daniel-VirtualBox:~$ sudo mkdir
[sudo] password for daniel:
mkdir: missing operand
Try 'mkdir --help' for more information.
daniel@daniel-VirtualBox:~$ sudo mkdir
mkdir: missing operand
Try 'mkdir --help' for more information.
daniel@daniel-VirtualBox:~$ mkdir documents
mkdir: cannot create directory 'documents': File exists
daniel@daniel-VirtualBox:~$ mkdir ramseylabwork
mkdir: cannot create directory 'ramseylabwork': File exists
daniel@daniel-VirtualBox:~$ cd documents
daniel@daniel-VirtualBox:~/documents$ ls
daniel@daniel-VirtualBox:~/documents$ cd Documents
bash: cd: Documents: No such file or directory
daniel@daniel-VirtualBox:~/documents$ cd ..
daniel@daniel-VirtualBox:~$ ls
181204_Promoter_Gpredict_antisense1.bed  examples.desktop  ramseylabwork
Desktop                                help              Templates
documents                             Music             Videos
Documents                             Pictures
Downloads                             Public
daniel@daniel-VirtualBox:~$ cd Documents/
daniel@daniel-VirtualBox:~/Documents$ ls
181115_LVSG4_Predictor_AntiSenseStrand.bed.txt  desktop.ini
181204_LVSG4_Predictor_SenseStrand.bed.txt      Promoters.bed
181204_Promoter_Gpredict_antisense1.bed         ramseylabwork
181204_Promoter_Gpredict_Sense1.bed
daniel@daniel-VirtualBox:~/Documents$

```

Correction from 18.12.05

```

bedtools intersect -wb -a Promoters.bed -b 181204_LVSG4_Predictor_SenseStrand.bed >
181204_Promoter_Gpredict_Sense1.bed

```

Preliminary Result: Could not open Exiting

I noticed the typo Strange instead of Strand for the file before the carrot. Success was obtained by correcting the typo and including the txt. file extension after the bed extension.

```

daniel@daniel-VirtualBox:~/Documents$ bedtools intersect -wb -a Promoters.bed -b
181204_LVSG4_Predictor_SenseStrange.bed> 181204_Promoter_Gpredict_Sense1.bed
Error: Unable to open file 181204_LVSG4_Predictor_SenseStrange.bed. Exiting.
daniel@daniel-VirtualBox:~/Documents$ bedtools intersect -wb -a Promoters.bed -b
181204_LVSG4_Predictor_SenseStrand.bed> 181204_Promoter_Gpredict_Sense1.bed

```

Error: Unable to open file 181204_LVSG4_Predictor_SenseStrand.bed. Exiting.

daniel@daniel-VirtualBox:~/Documents\$ bedtools intersect -wb -a Promoters.bed -b

181204_LVSG4_Predictor_SenseStrand.bed.txt 181204_Promoter_Gpredict_Sense1.bed

NC_007880	8212	8239	1	NC_007880	8212	8239	120
NC_007880	8259	8269	1	NC_007880	8259	8269	80
NC_007880	25495	25520	1	NC_007880	25495	25520	90
NC_007880	66340	66368	1	NC_007880	66340	66368	110
NC_007880	78869	78888	1	NC_007880	78869	78888	90
NC_007880	81216	81244	1	NC_007880	81216	81244	110
NC_007880	124548	124574	1	NC_007880	124548	124574	130
NC_007880	148724	148745	1	NC_007880	148724	148745	120
NC_007880	158523	158551	1	NC_007880	158523	158551	120
NC_007880	178059	178082	1	NC_007880	178059	178082	100
NC_007880	181637	181658	1	NC_007880	181637	181658	100
NC_007880	190574	190600	1	NC_007880	190574	190600	90
NC_007880	192525	192545	1	NC_007880	192525	192545	100
NC_007880	213661	213686	1	NC_007880	213661	213686	110
NC_007880	213945	213967	1	NC_007880	213945	213969	90
NC_007880	222996	223019	1	NC_007880	222996	223019	80
NC_007880	223338	223365	1	NC_007880	223338	223365	90
NC_007880	225226	225254	1	NC_007880	225226	225254	100
NC_007880	268117	268146	1	NC_007880	268117	268146	120
NC_007880	327029	327042	1	NC_007880	327029	327042	90
NC_007880	327770	327794	1	NC_007880	327770	327794	130
NC_007880	336032	336060	1	NC_007880	336032	336060	90
NC_007880	338822	338849	1	NC_007880	338822	338849	130
NC_007880	344375	344391	1	NC_007880	344375	344391	110
NC_007880	344055	344082	1	NC_007880	344055	344082	100
NC_007880	352392	352411	1	NC_007880	352392	352411	90
NC_007880	373707	373736	1	NC_007880	373707	373736	110
NC_007880	404281	404307	1	NC_007880	404281	404307	90
NC_007880	405453	405475	1	NC_007880	405453	405475	120
NC_007880	406551	406564	1	NC_007880	406551	406564	80
NC_007880	414280	414297	1	NC_007880	414280	414297	110
NC_007880	415569	415595	1	NC_007880	415569	415595	130
NC_007880	422922	422951	1	NC_007880	422922	422951	110
NC_007880	435559	435580	1	NC_007880	435559	435580	100
NC_007880	502373	502402	1	NC_007880	502373	502402	120
NC_007880	521242	521270	1	NC_007880	521242	521270	120
NC_007880	537122	537137	1	NC_007880	537114	537137	110
NC_007880	555874	555901	1	NC_007880	555874	555901	140
NC_007880	582517	582530	1	NC_007880	582517	582537	90
NC_007880	584676	584704	1	NC_007880	584676	584704	90
NC_007880	606167	606187	1	NC_007880	606167	606187	90
NC_007880	636615	636640	1	NC_007880	636615	636640	100
NC_007880	801864	801876	1	NC_007880	801864	801893	110
NC_007880	805457	805479	1	NC_007880	805457	805479	100
NC_007880	824484	824499	1	NC_007880	824484	824513	150
NC_007880	826501	826520	1	NC_007880	826501	826520	90
NC_007880	831838	831864	1	NC_007880	831838	831864	100
NC_007880	845448	845464	1	NC_007880	845442	845464	100
NC_007880	873852	873877	1	NC_007880	873852	873877	100
NC_007880	874836	874861	1	NC_007880	874836	874861	100
NC_007880	880022	880049	1	NC_007880	880022	880049	110
NC_007880	918322	918350	1	NC_007880	918322	918350	100
NC_007880	935709	935721	1	NC_007880	935709	935721	90
NC_007880	971339	971356	1	NC_007880	971339	971356	110
NC_007880	975531	975560	1	NC_007880	975531	975560	90

```

NC_007880 1015545 1015572 1 NC_007880 1015545 1015572 130
NC_007880 1028699 1028718 1 NC_007880 1028699 1028718 90
NC_007880 1029362 1029387 1 NC_007880 1029362 1029387 90
NC_007880 1055003 1055030 1 NC_007880 1055003 1055030 100
NC_007880 1086238 1086261 1 NC_007880 1086238 1086261 90
NC_007880 1119187 1119213 1 NC_007880 1119187 1119213 130
NC_007880 1130499 1130521 1 NC_007880 1130499 1130521 80
NC_007880 1144903 1144929 1 NC_007880 1144903 1144929 110
NC_007880 1172400 1172417 1 NC_007880 1172394 1172417 100
NC_007880 1184210 1184223 1 NC_007880 1184210 1184223 80
NC_007880 1184234 1184245 1 NC_007880 1184234 1184247 80
NC_007880 1199005 1199022 1 NC_007880 1199005 1199022 90
NC_007880 1201251 1201278 1 NC_007880 1201251 1201278 100
NC_007880 1229592 1229621 1 NC_007880 1229592 1229621 110
NC_007880 1232890 1232916 1 NC_007880 1232890 1232916 120
NC_007880 1232920 1232947 1 NC_007880 1232920 1232947 100
NC_007880 1282355 1282380 1 NC_007880 1282355 1282380 100
NC_007880 1293099 1293126 1 NC_007880 1293099 1293126 130
NC_007880 1323215 1323243 1 NC_007880 1323215 1323243 110
NC_007880 1338877 1338890 1 NC_007880 1338877 1338890 80
NC_007880 1338953 1338966 1 NC_007880 1338953 1338979 100
NC_007880 1340052 1340078 1 NC_007880 1340052 1340078 120
NC_007880 1349887 1349914 1 NC_007880 1349887 1349914 90
NC_007880 1420655 1420671 1 NC_007880 1420655 1420671 90
NC_007880 1437004 1437010 1 NC_007880 1437004 1437017 80
NC_007880 1463183 1463198 1 NC_007880 1463183 1463198 110
NC_007880 1490564 1490591 1 NC_007880 1490564 1490592 110
NC_007880 1522503 1522528 1 NC_007880 1522503 1522528 110
NC_007880 1572130 1572143 1 NC_007880 1572130 1572143 90
NC_007880 1648067 1648093 1 NC_007880 1648067 1648093 120
NC_007880 1648528 1648556 1 NC_007880 1648528 1648556 110
NC_007880 1683742 1683762 1 NC_007880 1683742 1683762 90
NC_007880 1685007 1685021 1 NC_007880 1685007 1685021 90
NC_007880 1686510 1686524 1 NC_007880 1686510 1686524 100
NC_007880 1701179 1701203 1 NC_007880 1701179 1701203 100
NC_007880 1746188 1746202 1 NC_007880 1746188 1746202 100
NC_007880 1847034 1847063 1 NC_007880 1847034 1847063 90
NC_007880 1861800 1861822 1 NC_007880 1861800 1861822 120
NC_007880 1881271 1881300 1 NC_007880 1881271 1881300 90
NC_007880 1881725 1881750 1 NC_007880 1881725 1881750 80
NC_007880 1884734 1884757 1 NC_007880 1884734 1884757 110
NC_007880 1894070 1894084 1 NC_007880 1894070 1894084 90

```

```

daniel@daniel-VirtualBox:~/Documents$ bedtools intersect -wb -a Promoters.bed -b
181204_LVSG4_Predictor_SenseStrand.bed.txt > 181204_Promoter_Gpredict_Sense1.bed

```

Viewing the header of the Promoter Gpredict_Sense file (first 10 lines)

```
cat 181204_Promoter_Gpredict_Sense1.bed | head
```

```

NC_007880 8212 8239 NC_007880 8212 8239 120
NC_007880 8259 8269 NC_007880 8259 8269 80
NC_007880 25495 25520 NC_007880 25495 25520 90
NC_007880 66340 66368 NC_007880 66340 66368 110
NC_007880 78869 78888 NC_007880 78869 78888 90

```

```

NC_007880 81216 81244 NC_007880 81216 81244 110
NC_007880 124548 124574 NC_007880 124548 124574 130
NC_007880 148724 148745 NC_007880 148724 148745 120
NC_007880 158523 158551 NC_007880 158523 158551 120
NC_007880 178059 178082 NC_007880 178059 178082 100

```

This was done to prepare to use the `wc -l` command

```

daniel@daniel-VirtualBox:~/Documents$ cat 181204_Promoter_Gpredict_Sense1.bed | wc -l
97

```

The significance of 97 is that there are purported to be 97 overlaps between the predicted Sense1 bed file and the reference Promoter bed file that contains experimentally determined promoters.

Original G4 predictions on LVS chromosome

```

daniel@daniel-VirtualBox:~/Documents$ cat 181204_LVSG4_Predictor_SenseStrand.bed.txt | wc -l
1293
daniel@daniel-VirtualBox:~/Documents$ cat 181115_LVSG4_Predictor_AntiSenseStrand.bed.txt | wc -l
1103

```

	Sense	Antisense	Total of both Strands
Overlap with Promoter file Prediction	97	169	266
# of Putative G4s	1,293	1,103	2,396
Percentage	0.07%	0.15%	0.11%

January 2019

Arbitrary PCR and Transposon Mapping

Monday, January 28, 2019

To Do:

- ~~1. Familiarize oneself with Arbitrary PCR~~
2. Create AR - PCR Protocol

Wednesday, January 30, 2019

To Do:

- ~~1. Create AR-PCR Protocol~~
- ~~2. Create AR-PCR Worksheet~~

Results and Data:

Delayed run of Arbitrary PCR since designed primers have yet to arrive.

PCR #1						
Component	Stock Concentration		Final Concentration		1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X		1X		5	25
Template	500ng		160ng/μL		3.13	indiv
dNTPs	10 mM		0.2mM		0.5	2.5

ddiH2O	-		-		14.245	71.225
--------	---	--	---	--	--------	--------

OneTaq Polymerase	0.25µL		1.25 units/50µL		0.125	0.625
KROL87 Arb1 or KROL 89 Arb 6	100µM		4.0µM		1	5
Internal Specific Primer: KROL90 Tn_Mar1	10µM		0.4µM		1	5
Total volume					25	

Run PCR samples through thermocycler program for reaction 1.

Initial

Denaturation: 5 min. 95°C

5-6X:

- Denature
30 sec.
94°C
- Anneal
30 sec.
30°C
- Extend
60 sec.
72°C

30X:

- Denature
30 sec.
94° C
- Anneal

30 sec.
45-55° C

- Extend
120 sec.
72°C

Final Extension:

- Extend 5
min.
72°C

PCR #2						
Component	Stock Concentration		Final Concentration		1 Rxn Volume (µL)	
PCR 1 Amplicon	–		–		0.75	3.75
dNTPs	10 mM		0.2mM		0.5	2.5
OneTaq Rxn Buffer	5X		1X		5	25
KROL88 Arb2	100µM		4.0µM		1	5
External Specific Primer: KROL91 Tn_Mar2	10µM		0.4µM		1	5
OneTaq Polymerase			1.25 units/50µL		0.125	0.625

ddiH ₂ O	–		–		16.625	83.125
Total volume					25	125

Run ArbPCR
reaction amplicons
through
Thermocycler
Program for
ARBPCR reaction
2.

PCR Rxn #2

30X:

- Denature:
30 sec.
94°C
- Anneal:
30 sec.
45-
55°C
- Extend:
60 sec.
72°C

Final Extension:

- 5 min.
72°C

February 2019

Monday, February 4, 2019

To Do:

- ~~1. Run Arbitrary PCR of 19/1/30~~

AP-PCR reactions 1 and 2 done with Arb1-Arb2 pairing subjected to 55°C annealing temperature during both reactions. No mixing of master mix components.

Wednesday, February 6, 2019



To Do:

1. Gel run Arbitrary Amplicons of 19.02.04

Ran a 1% agarose gel with 1X TAE buffer.

Well	Volume(uL)	Contents
1	10	Log ₂ DNA Ladder
2	30	PCR2 Tn3 Mutant
3	30	PCR2 Tn6 Mutant
4	30	PCR2 LVS

PCR samples volume samples were loaded as 25uL. 5uL of 6X orange G dye were added to each DNA template sample totaling a volume of 30uL. Not shown is the negative control



Results

Gel image does not depict any size bands less than 1kb, everything below 1kb has run off the gel. We are uncertain of the presence of any bands less than 1kb, and therefore have decided to retry PCR#2 from arbitrary PCR#1 stock. Currently the gel image shows bands present in the LVS lane as well as the PCR#2 Tn6 lane. LVS was our negative control of expected result of no bands. The presence of bands in the LVS lane indicates that amplification was arbitrary on the LVS gDNA. Retry scheduled for Friday February 8, 2019.

Friday, February 8, 2019

To Do:

1. Retry PCR#2 of 2/8/19
2. Retry Gel run of PCR #2 amplicons

Data:

Well	Volume(uL)	Contents
1	10	Log ₂ DNA Ladder
2	30	PCR2 Tn3 Mutant
3	30	PCR2 Tn6 Mutant

4	30	PCR2 LVS
---	----	----------

PCR #2

Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	
PCR 1 Amplicon	—	—	0.75	3.75
dNTPs	10 mM	0.2mM	0.5	2.5
OneTaq Rxn Buffer	5X	1X	5	25
KROL88 Arb2	100μM	4.0μM	1	5
External Specific Primer: KROL91 Tn_Mar2	10μM	0.4μM	1	5
OneTaq Polymerase		1.25 units/50μL	0.125	0.625
ddiH2O	—	—	16.625	83.125
Total volume			25	125

Run ArbPCR
reaction amplicons
through
Thermocycler
Program for
ARBPCR reaction
2.

PCR Rxn #2

30X:

- Denature:
30 sec.

94°C

- Anneal:
30 sec.
45-
55°C
- Extend:
60 sec.
72°C

Final Extension:

- 5 min.
72°C


This 2°Rxn of Arb1-Arb2 pairing was done at 55°C annealing temperature with no mixing of master mix components.

110mL of agarose in 1X TAE was poured into the gel casting tray and 6µL of SYBR Safe was added to the casting tray. A comb was used to help diffuse the SYBER Safe in the gel. I waited 25 minutes for the gel to solidify before loading samples and applying current. 10µL of 2log DNA ladder was added to lane 1 of the gel and 5µL of Orange G loading dye was added to each 25µL sample. Applied current was set for 90 mV until the gel front reached the midpoint of the gel. Successive runs from halfway down the gel were conducted at 90mV and 100-103mV with varied times. The gel was visualized using the gel dock reader and image lab a total of five times. Several images were taken but only four were saved. The saved images are the presented images below and are representative of different time points of the run.

Thoughts:

In comparing the gel images of this rerun through the experiment with the image of February 6, 2019 there is a striking difference in appearance and band presentation. Immediately clear are the bands in the No Template lane which suggest contamination. The image of February 6, 2019 shows bands with a lighting bolt appearance whereas the bands in the rerun gel are actual bands but much closer together, making size differential difficult but not impossible. When I was visualizing the gel using the gel dock reader and image lab I noticed that the gel was gliding along the surface of the dock platform when I was trying to position the gel within UV light focus. Gels are

like sponges especially when saturated with TAE buffer. Using Kim wipes I tried to soak up as much of the residual



TAE I could, to then take the picture.

Results:

Contamination in the No template lane voids all other results. Another rerun of this experiment is deemed appropriate.

Lanes from left to right are: 2log ladder, Tn3 mutant, Tn6 Mutant, LVS gDNA, and No Template




Image saved 08/02/19 1:22 PM


← The 2log DNA ladder in the image spans from 10kb to 0.5kb.

The Tn3 mutant lane shows two bands, one is of size 0.6kb and the other I think is 0.7kb.


Looking at the Tn6 mutant lane there is an intense band at the 0.8kb mark and a band matching up to 0.5 kb.

The LVS gDNA lane presents with one band that is less than 0.5kb. The 2log ladder in the image does not show any band less than 0.5kb but considering the proximity that the sole band in the LVS lane is to 0.5kb I consider it to be 0.4kb.


No template lane shows two intense bands. The bottom most lines up with the 0.4kb band found in the LVS lane. The intense band appears to be 0.8kb.



← This image is a duplicate of the image above with the brightness of the image raised and the fluorescent intensity of the bands highlighted through image lab. Image saved 08/02/19 1:24 PM



← This image is the result of running the gel further past the point of where the gel front reached the midpoint of the gel itself. Voltage current was set at 90V and time elapsed was 5 minutes. After 5 minutes I stopped the power source but did not immediately visualize the state of the gel. Image saved 08/02/19 1:38 PM



← The gel at this point was subjected to 100V of current for 7 minutes since the prior time it was visualized. Image saved 08/02/19 1:48 PM

Monday, February 11, 2019

To do:

1. Lab maintenance (Autoclave pipette boxes)
2. Lab Maintenance (Make Agar Media)

Wednesday, February 13, 2019

To do:

1. Retry PCR#2 Gel for Arbitrary PCR Transposon mutagenesis

PCR #1						
Component	Stock Concentration		Final Concentration		1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X		1X		5	25
Template	500ng		160ng/μL		3.13	indiv
dNTPs	10 mM		0.2mM		0.5	2.5
ddiH2O	-		-		14.245	71.225

OneTaq Polymerase	0.25μL		1.25 units/50μL		0.125	0.625
KROL87 Arb1 or KROL 89 Arb 6	100μM		4.0μM		1	5

Internal Specific Primer: KROL90 Tn_Mar1	10 μ M		0.4 μ M		1	5
Total volume					25	

Run PCR samples
through
thermocycler
program for
reaction 1.

Initial

Denaturation: 5

min. 95°C

5-6X:

- Denature
30 sec.
94°C
- Anneal
30 sec.
30°C
- Extend
60 sec.
72°C

30X:

- Denature
30 sec.
94° C
- Anneal
30 sec.
45-55° C
- Extend
120 sec.
72°C

Final Extension:

- Extend 5
min.
72°C

PCR #2						
Component	Stock Concentration		Final Concentration		1 Rxn Volume (μL)	
PCR 1 Amplicon	—		—		0.75	3.75
dNTPs	10 mM		0.2mM		0.5	2.5
OneTaq Rxn Buffer	5X		1X		5	25
KROL88 Arb2	100μM		4.0μM		1	5
External Specific Primer: KROL91 Tn_Mar2	10μM		0.4μM		1	5
OneTaq Polymerase			1.25 units/50μL		0.125	0.625
ddiH2O	—		—		16.625	83.125
Total volume					25	125

Run ArbPCR
reaction amplicons
through
Thermocycler

Program for
ARBPCR reaction
2.

PCR Rxn #2

30X:

- Denature:
30 sec.
94°C
- Anneal:
30 sec.
45-
55°C
- Extend:
60 sec.
72°C


Final Extension:

- 5 min.
72°C

AP-PCR done with Arb1 & Arb2 pairing at 55°C and with no mixing of master mix components to homogenize the solution.

Results:


Using the DNA Gel ladder from the Dutta Laboratory Hannah (HT) and I used the entirety of a 1% gel. HT loaded lanes 2-10 for her experiment. A lane was intentionally left blank to separate my samples from HT's samples. I loaded my samples in lanes 12-15. Lane 12 is transposon Tn3 mutant, Lane 13 is transposon Tn6 mutant, lane 14 is LVS gDNA, and finally lane 15 is No template. The gel shows a fluorescently strong band in the LVS lane and faint bands in all other lanes of interest.



Friday, February 15, 2019**To do:****1. Run Arbitrary PCR stage 1 with KROL 88 (Arb 6) primer****Results and Data:**

PCR #1				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X	1X	5	25
Template	500ng	160ng/μL	3.13	indiv
dNTPs	10 mM	0.2mM	0.5	2.5
ddiH2O	-	-	14.245	71.225
OneTaq Polymerase	0.25μL	1.25 units/50μL	0.125	0.625
KROL87 Arb1 (or KROL 88Arb 6)	100μM	4.0μM	1	5
Internal Specific Primer: KROL90 Tn_Mar1	10μM	0.4μM	1	5
Total volume			25	
PCR #2				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	

PCR 1 Amplicon	–	–	0.75	3.75
dNTPs	10 mM	0.2mM	0.5	2.5
OneTaq Rxn Buffer	5X	1X	5	25
KROL88 Arb2	100 μ M	4.0 μ M	1	5
External Specific Primer: KROL91 Tn_Mar2	10 μ M	0.4 μ M	1	5
OneTaq Polymerase		1.25 units/50 μ L	0.125	0.625
ddiH2O	–	–	16.625	83.125
Total volume			25	125



Bands in LVS lane and their intensity suggests many random primer hybridization events occurred. The Tn6 lane shows approximately similar size bands. These findings inspired the decision to image a PCR reaction 1 of samples 25.02.19. Lanes read left to right: DNA ladder, Tn3 mutant, Tn6 mutant, LVS gDNA, and No template.

Reaction occurred under 55°C annealing temperature and master with a master mix solution that were not thoroughly mixed well.


Monday, February 25, 2019

To do:

- 1. Run a 1% agarose gel with Arb6PCR reaction 1 samples.

Thoughts and Results

Wasteful errors during gel loading on my part caused poor quality gel with no information to be gleam.



Friday, March 1, 2019

To do:

- 1. Run Arb6 PCR reaction 1 and a 1% agarose gel from the generated products

March 2019

Monday, March 4, 2019

To do:

- 1.Run Arb6 PCR reaction 1 and a 1% agarose gel from the generated PCR products of Monday,

February 25, 2019

This AP-PCR trial was done with Arb6-Arb2 pairing at 55°C annealing temperature. No thorough mixing of the master mix components was done.

PCR #1						

Component	Stock Concentration		Final Concentration		1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X		1X		5	25
Template	500ng		160ng/μL		3.13	indiv
dNTPs	10 mM		0.2mM		0.5	2.5
ddiH2O	-		-		14.245	71.225

OneTaq Polymerase	0.25μL		1.25 units/50μL		0.125	0.625
KROL87 Arb1 or KROL 89 Arb 6	100μM		4.0μM		1	5
Internal Specific Primer: KROL90 Tn_Mar1	10μM		0.4μM		1	5
Total volume					25	

Run PCR samples through thermocycler program for reaction 1.

Initial

Denaturation: 5

min. 95°C

5-6X:

- Denature
30 sec.
94°C
- Anneal
30 sec.
30°C
- Extend
60 sec.
72°C

30X:

- Denature
30 sec.
94° C
- Anneal
30 sec.
45-55° C
- Extend
120 sec.
72°C

Final Extension:

- Extend 5
min.
72°C

PCR #2						
Component	Stock Concentration		Final Concentration		1 Rxn Volume (µL)	
PCR 1 Amplicon	—		—		0.75	3.75
dNTPs	10 mM		0.2mM		0.5	2.5

OneTaq Rxn Buffer	5X		1X		5	25
KROL88 Arb2	100 μ M		4.0 μ M		1	5
External Specific Primer: KROL91 Tn_Mar2	10 μ M		0.4 μ M		1	5
OneTaq Polymerase			1.25 units/50 μ L		0.125	0.625
ddiH2O	—		—		16.625	83.125
Total volume					25	125

Run ArbPCR
reaction amplicons
through
Thermocycler
Program for
ARBPCR reaction
2.

PCR Rxn #2

30X:


- Denature:
30 sec.
94°C
- Anneal:
30 sec.
45-
55°C
- Extend:
60 sec.

72°C

Final Extension:

- 5 min.
72°C

Results & Thoughts



Lanes Left to Right: DNA ladder, Tn3, Tn6, LVS gDNA, and No template. Seen here is AP-PCR primary reaction with Arb6-KROL90 pairing. Reaction done at 55°C and no mixing of master mix components.

Retry PCR reaction #1, this time describe every step of the preparation process and work process in as much detail as possible. Lack of band in the Tn3 mutant lane suggest that no PCR reaction occurred.

Wednesday, March 6, 2019

To do:

1. Rerun Arbitrary PCR reaction 1 with pKROL 89 (Arb6) primer

Protocol

Arbitrary PCR of embellished Transposon elements in mutant genomic or clonal DNA is a powerful technique that allows the researchers to map out where transposon elements have inserted into bacterial chromosome and allows answers to questions related to promoters' positions and phenotype profiling of mutants.

Create a PCR master mix for the first reaction

Samples

Sample	Tn3	Tn6	LVS gDNA	No Template
Concentration (ng/ul)	160	160	160	-

Arbitrary PCR Master Mix Reaction One Recipe

PCR #1				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X	1X	5	25
Template	500ng	160ng/μL	3.13	indiv
dNTPs	10 mM	0.2mM	0.5	2.5
ddiH2O	-	-	14.245	71.225
OneTaq Polymerase	0.25μL	1.25 units/50μL	0.125	0.625
KROL87 Arb1 (or KROL 88Arb 6)	100μM	4.0μM	1	5
Internal Specific Primer: KROL90 Tn_Mar1	10μM	0.4μM	1	5

Total volume			25	
--------------	--	--	----	--

Make up the master mix by aliquoting...

- i. 71.225uL of ddiH₂O
- ii. 25uL OneTaq Rxn Buffer
- iii. 2.5uL dNTPs
- iv. 5.0uL KROL 87 Arb1 or KROL 88 Arb6
- v. 5.0uL KROL 90 Tn_Mar1
- vi. 0.625uL OneTaq Polymerase

Once Master mix is prepared obtain four PCR strip tubes and label the four tubes as: Tn3, Tn6, LVS gDNA, and No Template.

Each PCR strip tube will contain 3.13uL of DNA sample and will total to 25uL. Load ~21.9uL of the master mix into each PCR strip tube and then aliquot 3.13uL of the DNA sample. Once the tubes contain master mix and DNA they are ready for the thermocycler.

Thermocycler Program

Initial Denaturation: 5 min. 95°C

5-6X:

- Denature 30 sec. 94°C
- Anneal 30 sec. 30°C
- Extend 60 sec. 72°C

30X:


- Denature 30 sec. 94° C
- Anneal 30 sec. 45-55° C
- Extend 120 sec. 72°C

Final Extension:

- Extend 5 min. 72°C

Once samples have been amplified run a gel electrophoresis using only 20uL of the amplicons. Melt 1% agarose gel in 1X TAE if the lab supplied with a flask. If not, make up 500ml of 1% Agarose in 1X TAE in a 1,000mL flask. Obtain an empty clean 1,000mL flask and pour 20ml of 25X TAE into it and fill it to 500mL with 480mL of ddiH₂O, add 5g of pure agarose powder, and finally add a stir bar. Heat and stir the flask until all the agarose powder is dissolved in solution. Pour ~100mL of lukewarm Agarose solution into the casting tray and add ~5uL of SYBERSafe, mix well with comb. Place the comb near the top of the gel to form wells and let gel sit for 25 minutes to solidify. Once the gel solidifies fill the two chambers until the gel is submerged with 1X TAE. Add 5uL of OrangeG DNA stain to 25uL samples. Add 10uL of DNA 1kbplus Ladder to well #1. The subsequent wells will have the stained DNA samples loaded onto them.

Thoughts and Results:



This primary reaction of Arb 6-KROL 90 pairing was subjected to 55°C annealing temperature.




With every trial of arbitrary PCR there is variability in the resultant gels of first PCR reaction. Nothing except the procedure is constant, sometimes bands appear and disappear in lanes and their sizes vary from one image to the next. In the image above there are two band smears in the Tn6 mutant lane and primer dimers in all lanes. How to validate this result is contingent on what is causing the variability of the gel images.

Friday, March 8, 2019

To do:

1. **Execute Arbitrary PCR with KROL88 Arb6 Primer and compare results to KMR's trail.**

In running this trial the samples present with consistency and symmetry. On reflection, since first undertaking this experiment endeavor I have not been spinning down the master mix or letting all components completely defrost. These results are based on careful prudent technique practices and the consistency between my samples and KMR is evident of no variability.



Results and Data:

Reactions during this trial ran under 55°C annealing temperature and master mix components and reagents were allowed to thaw and were thoroughly mixed well by vortex and sonication

Shared bands between Tn3, Tn6, and LVS are annealing products of same size what is noteworthy are the bands in lanes Tn3 and Tn6 that are above

Monday, March 18, 2019

To Do:**1. Execute ArbPCR1 with KROL87****Protocol**

Arbitrary PCR of embellished Transposon elements in mutant genomic or clonal DNA is a powerful technique that allows the researcher to map out where transposon elements have inserted into bacterial chromosome via RNA sequencing or ChIP Sequencing and allows answers to questions related to promoters' positions and phenotype profiling of mutants.

Create a PCR master mix for the first reaction

Samples

Sample	Tn3	Tn6	LVS gDNA	No Template
Concentration (ng/ul)	160	160	160	-

Arbitrary PCR Master Mix Reaction One Recipe

PCR #1				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X	1X	5	25
Template	500ng	160ng/μL	3.13	indiv
dNTPs	10 mM	0.2mM	0.5	2.5
ddiH2O	-	-	14.245	71.225
OneTaq Polymerase	0.25μL	1.25 units/50μL	0.125	0.625
KROL87 Arb1	100μM	4.0μM	1	5

Internal Specific Primer: KROL90 Tn_Mar1	10 μ M	0.4 μ M	1	5
Total volume			25	

Make up the master mix by aliquoting...

- vii. 71.225uL of ddiH₂O
- viii. 25uL OneTaq Rxn Buffer
- ix. 2.5uL dNTPs
- x. 5.0uL KROL 87 Arb1 or KROL 89 Arb6
- xi. 5.0uL KROL 90 Tn_Mar1
- xii. 0.625uL OneTaq Polymerase

Once Master mix is prepared obtain four PCR strip tubes and label the four tubes as: Tn3, Tn6, LVS gDNA, and No Template.

Each PCR strip tube will contain 3.13uL of DNA sample and will total to 25uL. Load ~21.9uL of the master mix into each PCR strip tube and then aliquot 3.13uL of the DNA sample. Once the tubes contain master mix and DNA they are ready for the thermocycler.

Thermocycler Program

Initial Denaturation: 5 min. 95°C

5-6X:

- Denature 30 sec. 94°C
- Anneal 30 sec. 30°C
- Extend 60 sec. 72°C

30X:

- Denature 30 sec. 94° C
- Anneal 30 sec. 45-55° C
- Extend 120 sec. 72°C

Final Extension:


- Extend 5 min. 72°C

Once samples have been amplified run a gel electrophoresis using only 20uL of the amplicons (Samples were run under 55°C annealing temperature). Melt 1% agarose gel in 1X TAE if the lab supplied with a flask. If not, make up 500ml of 1% Agarose in 1X TAE in a 1,000mL flask. Obtain an empty clean 1,000mL flask and pour 20ml of 25X TAE into it and fill it to 500mL with 480mL of ddiH₂O, add 5g of pure agarose powder, and finally add a stir bar. Heat and stir the flask until all the agarose powder is dissolved in solution. Pour ~100mL of lukewarm Agarose solution into the casting tray and add ~5uL of SYBERSafe, mix well with comb. Place the comb near the top of the gel to form wells and let gel sit for 25 minutes to solidify. Once the gel solidifies fill the two chambers until the gel is submerged with 1X TAE. Add 5uL of OrangeG

DNA stain to 25uL samples. Add 10uL of DNA 1kbplus Ladder to well #1. The subsequent wells will have the stained DNA samples loaded onto them.

1° Reaction occurred with 55°C annealing temperature and well thorough mixing of the master mix reagents and components.

Results:



Notes from KMR For 3/20/19

Template is PCR rxn 1 diluted 1:100 and 2nd PCR rxn annealing temperature @ 65 Celsius.

Wednesday, March 20, 2019

To Do: 1. Make a 1:100 of the Arbitrary PCR reaction 1 amplicons of 3/18/19 and run Arbitrary PCR reaction 2 with a melting temperature of 65°C

Thoughts/Reflection

1µL of Tn3, Tn6, LVS, and No template gDNA were used for a 1:100 to dilution leaving 9µL of each sample left over. Next, I loaded the diluted samples into the thermocycler and ran the following program...

30x 94°C 30"

45-55°C 30"

72°C 1'

72°C 5'

I did so under the assumption that once diluted the amplicons of reaction 1 that were made on 3/18/19, could at once be put through the thermocycler program for the second arbitrary PCR reaction since they are a solution of the PCR components: DNA polymerase, dNTPS, and primers. This is faulty thinking since Reaction 1 and Reaction 2 of the Arbitrary PCR protocol I have been following do not have the same set of annealing primers. Ultimately, my productivity on this day in laboratory amount to nothing more than a sham due to confusion and being swayed by an unsound assumption.

Friday, March 22, 2019

To Do: 1. Re-dilute 1:100 dilutions of Arbitrary PCR sample amplicons of 3/18/19 and suspend them into Arbitrary PCR#2 master mix aliquots, to run thermo cycler program for Arbitrary PCR reaction #2

This AP-PCR trial run occurred with 55°C annealing temperature and a thoroughly mixed master mix.

Protocol

Arbitrary PCR of embellished Transposon elements in mutant genomic or clonal DNA is a powerful technique that allows the researcher to map out where transposon elements have inserted into bacterial chromosome via RNA sequencing or ChIP Sequencing and allows answers to questions related to promoters' positions and phenotype profiling of mutants.

Samples

Sample	Tn3	Tn6	LVS gDNA	No Template
Concentration (ng/ul)	160	160	160	-


PCR #2				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (µL)	

PCR 1 Amplicon	–	–	0.75	3.75
dNTPs	10 mM	0.2mM	0.5	2.5
OneTaq Rxn Buffer	5X	1X	5	25
KROL88 Arb2	100µM	4.0µM	1	5
External Specific Primer: KROL91 Tn_Mar2	10µM	0.4µM	1	5
OneTaq Polymerase		1.25 units/50µL	0.125	0.625
ddiH2O	–	–	16.625	83.125
Total volume			25	125

Make up the master mix by aliquoting...

- i. 83.125uL of ddiH₂O
- ii. 25uL OneTaq Rxn Buffer
- iii. 2.5uL dNTPs
- iv. 5.0uL KROL 88 Arb2
- v. 5.0uL KROL 91 Tn_Mar2
- vi. 0.625uL OneTaq Polymerase

Results



Left to Right Lanes: DNA 1kb plus ladder, Tn3, Tn6, LVS gDNA, and Control(+). ArbPCR amplicon of KROL 89 Arb 6

Monday, March 25, 2019

To Do: Run arbitrary PCR protocol 1 with arbitrary primer KROL87 at 60°C annealing temperature



Lanes read Left to Right: DNA 1kb plus ladder, Tn3, Tn6, LVS gDNA, Control (+).

Results and Data:

The AP-PCR reaction master mix was thoroughly mixed well.

Arbitrary PCR Master Mix Reaction One Recipe

PCR #1				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X	1X	5	25
Template	500ng	160ng/μL	3.13	indiv

dNTPs	10 mM	0.2mM	0.5	2.5
ddiH ₂ O	-	-	14.245	71.225
OneTaq Polymerase	0.25µL	1.25 units/50µL	0.125	0.625
KROL87 Arb1	100µM	4.0µM	1	5
Internal Specific Primer: KROL90 Tn_Mar1	10µM	0.4µM	1	5
Total volume			25	

Make up the master mix by aliquoting...


- i. 71.225uL of ddiH₂O
- ii. 25uL OneTaq Rxn Buffer
- iii. 2.5uL dNTPs
- iv. 5.0uL KROL 87 Arb1 or KROL 89 Arb6
- v. 5.0uL KROL 90 Tn_Mar1
- vi. 0.625uL OneTaq Polymerase

Wednesday, April 3, 2019

To Do: Perform arbitrary PCR using KROL 87 and KROL 89 and run at 65°C annealing temperature under arb PCR thermocycler program #1

Results & Data:

The left portion of the gel image below is where the KROL 87 sample set was loaded (Lane 1 DNA 1kb plus ladder, Lane 2 Tn3 mutant, Lane 3: Tn6 mutant, Lane 4: LVS gDNA, and Lane 5: - control (H₂O)). The right side of the gel image is where KROL 89 sample set was loaded into with the same sample chronology as that of the sample set of lanes 1-5.



Monday, April 8, 2019

To Do: 1. Objective: Retry the ArbPCR 1° rxn of KROL 87 and KROL 89 at 65°C annealing temperature

Procedure:

Samples

Sample	Tn3	Tn6	LVS gDNA	No Template
Concentration (ng/ul)	160	160	160	-

Arbitrary PCR Master Mix Reaction One Recipe

PCR #1				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	Factor
OneTaq Rxn Buffer	5X	1X	5	25

Template	500ng	160ng/μL	3.13	indiv
dNTPs	10 mM	0.2mM	0.5	2.5
ddiH₂O	-	-	14.245	71.225
OneTaq Polymerase	0.25μL	1.25 units/50μL	0.125	0.625
KROL87 Arb1	100μM	4.0μM	1	5
Internal Specific Primer: KROL90 Tn_Mar1	10μM	0.4μM	1	5
Total volume			25	

Make up the master mix by aliquoting...

- vii. 71.225uL of ddiH₂O
- viii. 25uL OneTaq Rxn Buffer
- ix. 2.5uL dNTPs
- x. 5.0uL KROL 87 Arb1 or KROL 89 Arb6
- xi. 5.0uL KROL 90 Tn_Mar1
- xii. 0.625uL OneTaq Polymerase

Once Master mix is prepared obtain four PCR strip tubes and label the four tubes as: Tn3, Tn6, LVS gDNA, and No Template.

Each PCR strip tube will contain 3.13uL of DNA sample and will total to 25uL. Load ~21.9uL of the master mix into each PCR strip tube and then aliquot 3.13uL of the DNA sample. Once the tubes contain master mix and DNA they are ready for the thermocycler.

Thermocycler Program

Initial Denaturation: 5 min. 95°C

5-6X:

- Denature 30 sec. 94°C
- Anneal 30 sec. 30°C
- Extend 60 sec. 72°C


30X:

- Denature 30 sec. 94° C
- Anneal 30 sec. 45-55° C [65° C for increasing annealing specificity]
- Extend 120 sec. 72°C

Final Extension:

- Extend 5 min. 72°C

Thoughts & Results: Unlike Wednesday April 3, 2019 in which I concluded that I forgot to add the genomic DNA of transposons mutants Tn3, Tn6, LVS gDNA to the KROL 87 set of PCR tubes, this time I made certain to add 3.13 µL of gDNA to each PCR reaction tube. In hindsight I most likely did add the KROL 87 sample to their wells on 19.04.03. This is the second occurrence of no amplicons present in the KROL 87 set.



Wednesday, April 10, 2019

To Do: Perform Arb PCR 2° Reaction with 1:100 diluted amplicons of April 8, 2019 at 65° annealing temperature.

Procedure:

PCR #2				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (μL)	
PCR 1 Amplicon	–	–	0.75	3.75
dNTPs	10 mM	0.2mM	0.5	2.5
OneTaq Rxn Buffer	5X	1X	5	25
KROL88 Arb2	100μM	4.0μM	1	5
External Specific Primer: KROL91 Tn_Mar2	10μM	0.4μM	1	5
OneTaq Polymerase		1.25 units/50μL	0.125	0.625
ddiH2O	–	–	16.625	83.125
Total volume			25	125

Make up the master mix by aliquoting...

- i. 83.125uL of ddiH₂O
- ii. 25uL OneTaq Rxn Buffer
- iii. 2.5uL dNTPs
- iv. 5.0uL KROL 88 Arb2
- v. 5.0uL KROL 91 Tn_Mar2
- vi. 0.625uL OneTaq Polymerase


AP-PCR Thermo Cycler Program 2° Rxn

30x 94°C 30"

45-55°C 30" [65° C for increasing primer annealing specificity]

72°C 1'

72°C 5'



Friday, April 12, 2019

To Do: 1. Dilute AP-PCR from 4/8/19 1:10 and set up secondary PCR rxn upscaled to 50ul and run it under 65°C annealing temperature then PCR purify. Elute in 30ul of 0.1X EB buffer and run 5ul of each sample on gel.

Procedure:

PCR #2				
Component	Stock Concentration	Final Concentration	1 Rxn Volume (µL)	
PCR 1 Amplicon	–	–	0.75	3.75

dNTPs	10 mM	0.2mM	0.5	2.5
OneTaq Rxn Buffer	5X	1X	5	25
KROL88 Arb2	100µM	4.0µM	1	5
External Specific Primer: KROL91 Tn_Mar2	10µM	0.4µM	1	5
OneTaq Polymerase		1.25 units/50µL	0.125	0.625
ddiH2O	–	–	16.625	83.125
Total volume			25	125

Make up the master mix (50µL reaction) by aliquoting...

- i. 166.25uL of ddiH₂O
- ii. 50uL OneTaq Rxn Buffer
- iii. 5uL dNTPs
- iv. 10.0 uL KROL 88 Arb2
- v. 10.0uL KROL 91 Tn_Mar2
- vi. 1.25 uL OneTaq Polymerase

AP-PCR Thermo Cyclor Program 2° Rxn

30x 94°C 30"

45-55°C 30" [65° C for increasing primer annealing specificity]

72°C 1'

72°C 5'

QIAquick PCR Purification

Add 5 volumes buffer PB to 1 volume of the PCR reaction and mix. (50µL of PCR sample and 250 µL of PB buffer)

Place a QIAquick column into a 2ml collection/microcentrifuge tube.

Aliquot PB buffer-PCR sample into the spin column and centrifuge @ 13,000 rpm for 1 min. Discard flow-through.


Wash with 750µL of PE buffer to the QIAquick column and centrifuge @13,000 rpm for 1 min. Discard flow-through.

Centrifuge once more at 13,000 rpm for 1 min. to remove any residual wash buffer.

Elute DNA samples using 30µL of 0.1% EB buffer and let sit for 1 minute in the column. Centrifuge for 1 min. at 13,000 rpm.


Resolve the sample on 1% agarose gel and mix 1 µl Orange G dye to 5µL of sample.

Results:



Friday April 12, 2019

1:10 dilution of PCR secondary amplicons of 4/8/19



Wednesday April 10, 2019

1:100 dilution of PCR secondary amplicons 4/8/19

Tuesday, April 16, 2019

To Do: Quantify purified PCR products of 4/12/19 using the nanodrop and set up DNA sequencing reactions for Tn3 and Tn6 samples.

Results:

Sample ID	User ID	Date	Time	ng/ul	A260	A280	260/280
Tn3	Default	4/16/2019	1:36 PM	33.56	0.671	0.366	1.83
	0.312 0.008						
Tn6	Default	4/16/2019	1:37 PM	24.12	0.482	0.268	1.80
	0.255 0.026						

Monday, April 22, 2019

To Do:

1. Process sequencing results of 4/19/19 through BLAST and record genetic loci of Tn3 and Tn6 samples. Next figure out where transposon insertion and flanking occurred. Annotate the insertion and flanking on the LVS genome file of Snappgene.
2. Create KROL primers: For Tn3 create primers that anneal 100 bp upstream and downstream from Transposon. For Tn6 create primers that anneal 150 bp upstream and downstream from the transposon.

Results:

Genbank information for Tn3 sample

Francisella tularensis subsp. holarctica LVS complete genome

Sequence ID: [AM233362.1](#) Length: 1895994 Number of Matches: 1

Range 1: 698190 to 698244 [GenBankGraphics](#) Next Match Previous Match

Alignment statistics for match #1

Score Expect Identities Gaps Strand

102 bits(55) 2e-23 55/55(100%) 0/55(0%) Plus/Plus

Query	1	AAAGTACCACTAGTTTTAGATGTAATATATTTTCTATGTCGGGATTATCTTCAT	55
Sbjct	698190	AAAGTACCACTAGTTTTAGATGTAATATATTTTCTATGTCGGGATTATCTTCAT	698244

Genbank Information for Tn6 sample

Francisella tularensis subsp. holarctica LVS complete genome

Sequence ID: [AM233362.1](#) Length: 1895994 Number of Matches: 1Range 1: 721139 to 721198 [GenBankGraphics](#) Next Match Previous Match

Alignment statistics for match #1

Score	Expect	Identities	Gaps	Strand
111 bits(60)	3e-26	60/60(100%)	0/60(0%)	Plus/Plus
Query 1		CACATCAAAAAGTTGAGATATCAGGTATACCATATTACCGATAATCACAGCTTTCAAAGC		60
Sbjct 721139		CACATCAAAAAGTTGAGATATCAGGTATACCATATTACCGATAATCACAGCTTTCAAAGC		721198

Wednesday, April 24, 2019**To Do: 1. Using the newly created primers of 4/22/19 validate findings of 4/22/19 with PCR****Procedure**

Total reaction volume	20
Total number of reactions	4

				Factor
Component	Stock concentration	Final concentration	1 rxn volume	5.5
ddiH2O			4	22
KOD buffer	2x	1x	10	55
dNTPs	2 mM	0.4 mM	4	22
oligo F	10 uM	0.3 uM	0.6	3.3
oligo R	10 uM	0.3 uM	0.6	3.3
template	100 ng/ul	2 ng/ul	0.4	
KOD	1 U/ul	0.02 U/ul	0.4	2.2
Total volume			20	110

Make Two Master Mixes

1. MM#1 KROL 160 + KROL 161
2. MM#2 KROL 162 + KROL 163

KROL 160 + KROL 161

Tn size = 4616

WT KROL 160 and KROL 161 = 221 bp

WT product bands = (4616 + 221) = 4837 bp

KROL 162 + KROL 163

Tn size = 4616 bp

WT KROL 162 and KROL 163 = 358 bp

WT product bands = (4616 + 358) = 4974 bp

PCR Program

Heat at 94 degrees for 2 minutes,

94 degrees C for 20 seconds

50 degrees C for 30 seconds

68 degrees C for 1 minute (KOD polymerase functions properly at 68 degrees C; Taq polymerase is different temp)


Go back to step 2

Repeat 32x

68 degrees C for 5 minutes

12 degrees C for infinity

Results: Failing to set the extension time for 5 minutes during the multiple cycles of PCR gave low abundance of amplicon products. KOD polymerase requires time for nucleotide polymerization. Running the gel for close to an hour allowed a portion of the gel front to run off the gel.



Friday, April 24, 2019

To Do:

1. Run Tn3, Tn6, LVS gDNA, and H₂O PCR reactions with KROL 160 and KROL 161 primers.

Procedure

Total reaction volume	20
Total number of reactions	4

				Factor
Component	Stock concentration	Final concentration	1 rxn volume	5.5
ddiH ₂ O			4	22
KOD buffer	2x	1x	10	55
dNTPs	2 mM	0.4 mM	4	22
oligo F	10 uM	0.3 uM	0.6	3.3
oligo R	10 uM	0.3 uM	0.6	3.3
template	100 ng/ul	2 ng/ul	0.4	
KOD	1 U/ul	0.02 U/ul	0.4	2.2

Total volume	20	110
--------------	----	-----

Make One Master Mix using KROL 160 and KROL161

1. **MM#1 KROL 160 + KROL 161**
2. **MM#2 KROL 162 + KROL 163**

KROL 160 + KROL 161

Tn size = 1441 bp

WT KROL 160 and KROL 161 = 221 bp

WT product bands = (1441 + 221) = 1662 bp

KROL 162 + KROL 163

Tn size = 1441 bp

WT KROL 162 and KROL 163 = 358 bp

WT product bands= (1441+ 358) = 1799 bp

PCR Program

Heat at 94 degrees for 2 minutes,

94 degrees C for 20 seconds


50 degrees C for 30 seconds

68 degrees C for 5 minutes (KOD polymerase functions properly at 68 degrees C extending elongation step to 5 minutes grants KOD polymerase a greater turnover number) Go back to step 2

Repeat 32x

68 degrees C for 5 minutes

12 degrees C for infinity

Results:

Originally on April 24th, 2019 the estimated size value of the transposon inserts was 4616bp, it wasn't until KMR realized the mishap of using the wrong reference plasmid map for results verification and results were evaluated with another reference plasmid map estimating inserts to be around 1.7 kb.

Bibliography

Ramsey, K. M. and Dove, S. L. (2016) 'A response regulator promotes *Francisella tularensis* intramacrophage growth by repressing an anti-virulence factor', *Molecular Microbiology*. doi: 10.1111/mmi.13418.