

Table of Contents

January 2020	2
Future To-Do	4
Bibliography	5

January 2020**Friday, January 10, 2020****To Do:**

1. Run gel for antibody tests
2. Start O/N transfer for antibody tests

Results and Data:

Run 2x 4-12% Bis-Tris gels using 1x MOPS buffer, with the following loading / plan:

1	LVS (1, 12.14.19)	Tul4	1	LVS (1, 12.14.19)	PdpD
2	Δ rpsU2 (3, 12.14.19)		2	Δ rpsU2 (3, 12.14.19)	
3	Magic Mark ladder		3	Magic Mark ladder	
4	LVS (1, 12.14.19)	IglA	4	LVS (1, 12.14.19)	PdpC
5	Δ rpsU2 (3, 12.14.19)		5	Δ rpsU2 (3, 12.14.19)	
6	Magic Mark ladder		6	Magic Mark ladder	
7	LVS (1, 12.14.19)	IglB	7	LVS (1, 12.14.19)	PdpB
8	Δ rpsU2 (3, 12.14.19)		8	Δ rpsU2 (3, 12.14.19)	
9	Magic Mark ladder		9	Magic Mark ladder	
10	LVS (1, 12.14.19)	IglC	10	LVS (1, 12.14.19)	PdpA
11	Δ rpsU2 (3, 12.14.19)		11	Δ rpsU2 (3, 12.14.19)	
12	Magic Mark ladder		12	Magic Mark ladder	
13	LVS (1, 12.14.19)	IglD	13	LVS (1, 12.14.19)	
14	Δ rpsU2 (3, 12.14.19)		14	Δ rpsU2 (3, 12.14.19)	
15	Magic Mark ladder		15	Magic Mark ladder	

Ran 10 μ L of each lysate (from pellets grown in sBHIc for T6SS assay by HT) and 1 μ L ladder

Used new gel running chambers. Ran gels at 150V for ~1 hr.

Transfer using new blot modules. Put at 4°C O/N, running at 7V (start ~5:20pm). New module does not fit in ice bucket to keep cool.

Saturday, January 11, 2020**To Do:**

1. Stop transfer, put membranes in blocking buffer at 4°C

Results and Data:

After ~ 17 hrs of transfer (7V, 4°C), stopped transfer. Cut membranes to pieces with 3 well each (hopefully!) and incubate in blocking buffer in cold room, rocking, until Monday.

Wednesday, January 15, 2020**To Do:**

1. Perform WB

Results and Data:

Blot	Protein	kDa	Rec dilution	Planned dilution	Volume in 10 mL
1	Tul4	17	1:5000	1:20000	0.5
2	IglA	20.9	1:2000	1:2000	5
3	IglB	57.9	1:1000	1:1000	10
4	IglC	22.1	1:1000	1:1000	10
5	IglD	46.4	1:1000	1:1000	10
6	PdpC	155.9	1:2500	1:1000	10
7	PdpB	127.5	1:1000	1:1000	10
8	PdpA	95.3	1:1000	1:500	20
9	PdpB	127.5	1:1000	1:500	20
10	PdpA	95.3	1:1000	1:200	50

Don't bother with PdpD antibody- would be 140kDa, but is a pseudogene in LVS.

Incubate in 10 mL 1° for 1 hr

Wash 4 x 10' (approx.- sometimes a bit more)

Re-block in 5 mL blocking buffer for 20' (+ probably 10-15' more)

Add 2° antibody directly to 5 mL blocking buffer- 1 uL of goat anti-mouse or goat anti-rabbit, as indicated below. Incubate 1 hr.


Blot	Protein	Antibody type
1	Tul4	mouse monoclonal
2	IgIA	rabbit polyclonal
3	IgIB	mouse monoclonal
4	IgIC	mouse monoclonal
5	IgID	mouse monoclonal
6	PdpC	rabbit polyclonal
7	PdpB	mouse monoclonal
8	PdpA	mouse monoclonal
9	PdpB	mouse monoclonal
10	PdpA	mouse monoclonal

Wash 4 x 10'

Develop ~3' with 5 mL developer mix each

Initial general impressions:

Too much anti-rabbit secondary. Both blot 2 and 6 have LOTS of background
Anti-PdpA has a background band (just as here, in the CoA):



Actually, I thought the background band on my blot was above the PdpA protein...

May need to look more carefully at the ladder. Should be at 93 kDa, the background on this blot is closer to 70 kDa.

Can't really see PdpC. Maybe if there was less background? The CoA blot isn't particularly convincing.

IgIB, IgIC, IgID antibodies look great!

IgIA looks pretty good with a VERY brief exposure. Will probably be great with less background.

Tul4 / LpnA continues to look excellent.

WT vs Δ rpsU2 impressions

Much less PdpA!

Much less PdpB!

Can't tell about PdpC

May actually be slightly less IgIB, C, D – although this may be a loading issue. No control, and LpnA/Tul4 blot is cut off, so can't tell relative amounts.

Aliquot plan

Tul4 / LpnA : 3 uL

IgIA: 10 uL

IgIB: 20 uL

IgIC: 20 uL

IgID: 20 uL


PdpC: 50 uL

PdpB: 20 uL

PdpA: 40 uL

Made 5 aliquots of each (only 4 of Tul4- there is only 100 uL total and I'm sending 15 uL to Mougous lab)

Scans of the blots:



Note that for PdpA, there is an extra band running around the size of PdpB, which is 127.5 kDa. It makes sense that the slightly lower band (closer to the 100 kDa band) would be PdpA, which is 95.3 kDa.

Future To-Do

Move 1° LVS pKR10-1 into strain box

Bibliography