

## Table of Contents

November 2021 .....	2
Lab goals.....	2
Future To-Do .....	6

## November 2021

### Lab goals

Paper: PmrA / PriM

Re-do macrophage assay with  $\Delta$ pmrA strains

Isolate gDNA from  $\Delta$ pmrA experiment in (April?) and check by qPCR

Paper:  $\Delta$ rpsU2

Re-do macrophage assay with  $\Delta$ rpsU2 and complements

### Sunday, November 7, 2021

#### To Do:

1. Prep fresh DMEM-F
2. Thaw macrophage (P5)

#### Results and Data:

Supplemented more DMEM with 10% FBS.

Thawed 1 vial of J774A.1 P5 cells, seeded into a 100 mm TC dish.

### Monday, November 8, 2021

#### To Do:

1. Feed macrophage
2. Pour plates

#### Results and Data:

Replaced media (removed 10 mL, washed with 5 mL, added back 10 mL DMEM-F). Cells look happy under microscope.

Dan and Aisling helped me autoclave 4 flasks of CHA. I poured plates for the macrophage assay for Thursday

**Wednesday, November 10, 2021****To Do:**

1. ~~Streak cells for mac expt~~
2. ~~Seed macrophage~~

**Results and Data:**

Washed cells with 5 mL DMEMF, added back 8 mL and scraped up cells. Pelleted at 1,000xg for 5 mins in 15 mL conical tube- good pellet! Resuspended in 7mL DMEMF, plenty of cells (want at least 12 mL at  $1.25 \times 10^5$ ).

1st measurement			
78	63	56	62
Average			64.75
Undiluted			129.5
Density			$1.30 \times 10^6$

Dilute a bit high just in case (better to have more than less)- use 1.5 mL + 13.5 mL DMEMF. Dispense using multichannel into 2x 96 well plates and check concentration:

2nd measurement			
5	8	7	3
8	5	6	8
5	10	6	5
Average			6.333333333
Undiluted			12.66666667
Density			$1.27 \times 10^5$

Perfect!

Use all remaining diluted cells and some of the resuspended cells (to make a volume of 10mL) and put into another TC dish as P6 cells. Put 96 well plates and dish into 37°C incubator with 5% CO<sub>2</sub>. Done around 5pm.

Final calculations in detail:

**Macrophage Calculations**

Cells per well	$2.50 \times 10^4$
Volume to plate (mL)	0.2
Density needed (cells/mL)	$1.25 \times 10^5$
Total volume needed (mL)	15
Measured cells per ml	$1.30 \times 10^6$
Volume stock needed (mL)	1.5
Volume media for dilution	15.0
Measured cells per ml, seeded	$1.27 \times 10^5$
Measured cells per well	$2.53 \times 10^4$

**Thursday, November 11, 2021**

**To Do:**

- 1. ~~Infect macro~~
- 2. ~~Plate inoculum~~
- 3. ~~Wash macro~~
- 4. ~~T=2~~

**Results and Data:**

Plan for inoculum:

Bacterial Calculations	Actual
MOI	5
Macrophage cells per well	2.53E+04
Volume bacteria to add (mL)	0.05
Bacterial density needed (cells/mL)	2.53E+06
Cells/mL per OD600	5.81E+09
OD needed for given density	0.00044
Resuspend to	0.050
Final MOI 5, dilute 1:100	0.00050

Resuspend bacteria, check ODs, dilute and re-check, perform final dilution for a final volume of 4 mL in 3-chamber multi reservoirs.

Number	Strain	Resuspend cells to (OD600)	For final vol 1 mL at 0.05		For final vol 4 mL at 0.0005	
			Cells (uL)	OD600	Cells (uL)	Volume media (uL)
1	LVS pF	2.30	21.7	0.049	40.0	3960
2	$\Delta$ rpsU2 pF	1.23	40.7	0.049	40.0	3960
3	$\Delta$ rpsU2 pF-rpsU1	2.00	25.0	0.051	40.0	3960
4	$\Delta$ rpsU2 pF-rpsU2	2.48	20.2	0.050	40.0	3960
5	$\Delta$ rpsU2 pF-rpsU3	2.17	23.0	0.054	39.0	3961
6	$\Delta$ pigR	3.19	15.7	0.052	40.0	3960

Add 50 ul each inoculum to macrophage as below, 2x, put macrophage infection plates back in incubator. **Infected at 11:40am.**

	1	2	3	4	5	6	7	8	9	10	11	12
A	LVS pF	LVS pF	LVS pF		$\Delta$ rpsU pF	$\Delta$ rpsU pF	$\Delta$ rpsU pF		$\Delta$ rpsU pF-rpsU1	$\Delta$ rpsU pF-rpsU1	$\Delta$ rpsU pF-rpsU1	
B												
C												
D	$\Delta$ rpsU pF-rpsU2	$\Delta$ rpsU pF-rpsU2	$\Delta$ rpsU pF-rpsU2		$\Delta$ rpsU pF-rpsU3	$\Delta$ rpsU pF-rpsU3	$\Delta$ rpsU pF-rpsU3		$\Delta$ pigR	$\Delta$ pigR	$\Delta$ pigR	
E												
F												
G												
H												

Serially dilute inoculums 10x in PBS in 96-well plate, plate  $10^0 - 10^4$  on square CHA-Kan (#1-5) or square CHA (#6). Put plates at 37°C.

at 1:40pm, remove media from wells. Wash 2x with 1x PBS (280 uL) and replace media with 200 ul DMEMF with 10 ug/mL gentamycin.

$$10 \text{ ug/mL} * 12,000 \text{ ul} / 50,000 \text{ ug/mL} = 2.4 \text{ ul}$$

At 3:40, remove 1 plate from incubator. Remove media from wells. Wash 2x with 1x PBS (250 uL) and replace media with 200 ul 1% saponin in 1xPBS.

**Friday, November 12, 2021****To Do:**

1. T=24

**Results and Data:**

Words go here.

**Saturday, November 13, 2021****To Do:**

1. Check / pull plates

**Results and Data:**

Can count plates with inoculums from 1, 3, 4, and 6 – pull, count, and leave on bench

**Sunday, November 14, 2021****To Do:**

1. Check / pull plates

**Results and Data:**

Can count plates with inoculums from 2 and 5 – pull, count, and leave on bench  
Pull and leave on bench the following T=2 plates: 1, 4, 6

**Monday, November 15, 2021****To Do:**

1. T=24

**Results and Data:**

Pull and count the following T=2 plates: 2, 3  
Pull and count the following T=24 plates: 1, 2 (2-1 and 2-2; 2-3 needs another day), 3, 4

## Future To-Do

Move 1° LVS pKR10-1 into strain box

## Bibliography

Suh, Moo-Jin et al. "Extending ribosomal protein identifications to unsequenced bacterial strains using matrix-assisted laser desorption/ionization mass spectrometry." *Proteomics* vol. 5,18 (2005): 4818-31. doi:10.1002/pmic.200402111