

P3C Plasmid 3 -- BSL2

P3C Plasmid 3 -- BSL2

**Order Confirmation**

Kathryn Ramsey, kramsey@uri.edu

Order type: **Plasmid**  
**BSL2****3** samples, Order Code: **P3C**

Thank you for choosing plasmidsaurus!

Please make sure your samples EXACTLY MATCH the info entered in your order details. Discrepancies may delay or prevent sample processing.

Samples should ONLY be sent in **transparent, 96-well, PCR strip tubes** (or single PCR tubes if shipping single samples) or **96-well PCR plates**. Other type of tubes (1.5 mL, 0.5 mL, 384-well format, opaque, etc.) are not compatible with our rapid workflow.


Optionally, on the back of the paper, draw a dinosaur.

If our kids like it, you might win a prize. Art becomes our sole property to exploit as we wish.


**Shipping Instructions****STEP 1: Label your sample with your 3-character order code**

Each order will be assigned a **3-character order code** after your order has been confirmed through your plasmidsaurus account.

- Write your **order code** on the 1st tube of each strip.
- Write the **strip number** on the 2nd tube of each strip.



- Please label your tubes with the **order code** and **strip number** only. Do not include any other information, such as sample names.
- If you're sending a single sample in a single PCR tube, label with the 3-character order code only.
- Array samples in COLUMNS, NOT ROWS in 96-well plates, as listed in your order details.

**STEP 2: Package your sample****If sending PCR strip tubes or single PCR tubes:**

- Ensure caps fit snugly and close them firmly.
- Wrap caps in one light layer of parafilm to keep the caps closed.
- Place the wrapped tubes into a small protective vessel (such as a Falcon tube or small box).

**If sending PCR plates:**

- Plates must be sealed extremely thoroughly to survive shipping without leaking.
- We recommend using a heat seal (rather than adhesive) or strip caps on top of the plate.
- Wrap each plate individually in bubble wrap or other packaging material to prevent cracking.

Room temperature shipping is fine, dry ice and cold packs are not necessary.

**STEP 3: Ship your sample along with order information**

Make sure to include a copy of the **order confirmation** received upon submitting your order with your shipment.

You can choose to send us your order directly **OR** use one of our dropboxes:

- **Send order directly to us**

Samples should be sent via express courier (UPS, FedEx, DHL, etc.) to the following address:

**plasmidsaurus**  
1850 Millrace Drive, Suite 200  
Eugene, OR 97403  
Phone #: 541-514-0871  
Monday - Friday: 9 AM - 5 PM (PST)

- **Use our dropbox**

Simply place your order in our dropbox and enjoy free overnight shipping! No shipping envelope required, just put your samples and order info in a bag or Falcon tube.

P3C Plasmid 3 -- BSL2

P3C Plasmid 3 -- BSL2