

doi:10.1101/pdb.rec8297

Cold Spring Harb. Protoc. 2006. 2006: pdb.rec8297-

Recipe

SSC

For a 20X solution: Dissolve 175.3 g of NaCl and 88.2 g of sodium citrate in 800 mL of H₂O. Adjust the pH to 7.0 with a few drops of a 14 N solution of HCl. Adjust the volume to 1 L with H₂O. Dispense into aliquots. Sterilize by autoclaving. The final concentrations of the ingredients are 3.0 M NaCl and 0.3 M sodium citrate.